Homogeneização e retirada de bolhas de revestimento
Homogenization and withdrawal of covering bubblesh
Homogeneización y retirada de burbujas de revestimentos

turbomix

EDG
EQUIPAMENTOS

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turbomix

Homogenization and withdrawal of covering bubblesh

Date of manufacture

		200
--	--	-----

Serial Number

--

1 Characteristics

- Integrated building which incorporates all the elements in a single practical and strong cabinet.
- It can be used on table base or fixed on the wall.
- A single motor which activates both the vacuum pump and the mixer.
- Specially built high capacity vacuum pump does not need lubrication, it quickly removes the bubbles of the covering.
- Primary and secondary filters which eliminate the sludge that can harm the vacuum pump operation.
- Check valve which keeps the vacuum even after the spatulating process.
- 400RPM mixer which efficiently homogenizes the covering.
- Exclusive tank bracket which keeps the tank in its position during spatulation.
- Timer which determines the process time.
- Vacuumeter which indicates the vacuum level.
- 2 complete tanks, a table base, a wall bracket, bolts, adapter sleeves and instruction manual go with the equipment.
- 12 months warranty for manufacturing defects.

2 Reception

When receipting your TURBO MIX check its general condition. Check the components that go with the equipment:

- Main unit.
- Two tanks with lids.
- Wall fixing plate.
- Two screws with wall anchors
- Table base.
- Instruction manual.
- Warranty certificate.
- Consumer questionnaire.

Warning: Contact the company in case of any irregularity or lack of any of the components. Do not forget to fill in the consumer questionnaire and send it to the company.

3 Installing

There are two ways for installing TURBO MIX:

- 1- Fixed to a wall: Use the fixing plate supplied with the equipment and, with screws and anchors, fix the plate to the wall at an appropriate position.
- 2- Fixed to the base that goes with the equipment.

In both cases, mount the plate protuberances either on the base or on the equipment holes.

- Install your equipment at an appropriate place for its utilization.
- Check if the outlet voltage in which TURBO MIX will be plugged in is the same as the one indicated on the switch at the bottom of the equipment. TURBO MIX is shipped with 220v already selected.

For your own safety, do not forget to connect the grounding to the plug third pin (round pin), for doing so NEVER use the electrical network neutral, if you have any doubts contact an electrician.

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4 Testing

Once your TURBO MIX is installed it is ready to be tested.

- Turn the TEMPO button from right to left, as indicated by the arrow, and set a 30 seconds time.
- Cover the vacuum pen hole with your thumb.
- Turn the main switch on.
- Your TURBO MIX will start to work and the vacuum level will indicate the vacuum level (with a minimum of 24.5 inches). Once the time is over the equipment will shut down and the vacuum level will decrease slowly.

5 Use

- Set the spatulation time as indicated by the manufacturer of the covering to be used.
- Prepare the liquid covering pre-mixture using one of the tanks that goes with the equipment.
- Clean out all the covering remains that eventually stick on the tank edge.
- Adapt the tank lid to it.
- Mount the vacuum pen on its base on the tank lid.
- Lower the tank base and mount it.
- Turn the equipment on.
- When the time is over, lower the base, remove the tank and disconnect the vacuum pen.

6 Maintenance

Your TURBO MIX practically does not need maintenance, but some special attention should be given to its tanks and lids which should be washed with running water just after usage to avoid covering accumulation which, once hardened, is very difficult to be removed. The accumulation of covering generates the entrance of air in the tank or the blocking of the beater.

Clean the filters according to the instructions printed on the equipment from time to time. Do not turn the covering under vacuum upside down, this will cause the instantaneous clogging of the primary filter.

Your TURBO MIX vacuum pump was designed using technology which allows high performance without lubricants.

Therefore, it is not needed nor recommended the equipment lubrication.

EDG Equipamentos will always be ready to attend you in case your TURBO MIX needs settings or repair. However, even if the equipment is under warranty, before asking for technical assistance check if the problem doesn't have a simple solution. This can avoid waste of time and unnecessary removal of the equipment.

PROBLEMS	SOLUTIONS
The equipment does not work.	<p>Check if there is tension in the electrical network outlet.</p> <p>Check if the plug is well connected.</p> <p>Check if the commuter switch is in the correct position.</p> <p>Check if the timer button indicates "0" (zero).</p>
When the vacuum pen hole is covered with the thumb, the vacuum level does not reach the minimum level of 24.5 in/hg	<p>Check the vacuum pen tightness and sealing ring.</p> <p>Check the secondary filter body tightness and sealing ring.</p>
The vacuum level does not reach the minimum value during spatulation.	<p>Check if the vacuum pen is correctly adapted to its base.</p> <p>Check if there is not covering remains on the tank lid sealing ring and on its base.</p> <p>Check if there are not crushed parts on the tank which impede a good sealing of the lid.</p>
The vacuum level indicates an ideal vacuum level but there are bubbles in the covering.	<p>The primary or secondary filters are clogged obstructing the air passage from the tank to the pump, clean the primary and secondary filters.</p> <p>The vacuum pen base hole on the tank lid is clogged, clean it.</p>

SPECIFICATIONS:

110 and 220 volts (commuter switch), 180 watts, 50/60hz.
 Minimum vacuum 24.5 in/hg (inches of mercury).
 Timer 0 – 70 sec.
 Spatulator rotation 400 rpm (60Hz).
 Complete net weight 18 kg.
 Gross weight 22 kg.

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