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INSTRUCTION MANUAL

press²
alumini

Furnace for conventional and pressed ceramic.

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		200
--	--	-----

Número de Série

--

1 Recepting

The equipment is packed separately from the vacuum pump, check if you have two boxes if you have acquired the furnace and the pump. Check the general conditions of the package, if any damage is noticed contact the company immediately which will provide guidance on how to proceed. We remind you that the equipment is shipped at the purchaser's account and risk and it is insured by the conveyer.

1.1 The Sinter Press package should include:

- A) A furnace settled in injected foam.
- B) A spool for conventional ceramic Pic4-Item-14.
- C) A spool for pressed ceramic. Pic-5-item-17.
- D) A rigid mat -I – Pic-4-Item-15
- E) A soft mat.
- F) A refractory platform for pressed ceramic. Pic-5-Item-18.
- G) A hollow refractory platform for ceramic

- H) A set of isothermal pins. Pic-4-item-16
- I) Three 1.5 meters N. A. T. hoses
- J) Metallic clamps. (6 pieces).
- K) Instruction manual.
- L) Warranty deed.

1.2 The vacuum pump package should include:

- A) The vacuum pump
- B) Rubber, floor insulating caps for fixation (4 pieces).
- C) Labels for connections identification.

2 Main characteristics

- Quartz muffle with thermal insulation ceramic mat molded through low thermal mass vacuum, which insures perfect temperature homogeneity.
 - 79 programs: 78 for conventional and pressed ceramic and one standard recipe that comes from the factory.
 - Exclusive tray movable system to facilitate the positioning of works done on the platform.
 - Management by state-of-the-art microprocessing.
 - Easy operation and visualization of the burn parameters.
 - Ambient working temperatures at 1,100°C.
 - Maximum temperature protection limiter at 1,200°C
 - Linear heating speed from 0.1 to 70.0°C/min.
 - Thyristored power control
 - Burn time for conventional/pressed ceramic from 0 to 99 minutes.
 - Decreasing time visualization for all processes.
 - Vacuum time from 0 to 99 minutes with decreasing time visualization.
 - Adjustable and independent ascent and descent elevator times.
 - Safety systems protect the equipment from operation mistakes or defects.
 - Flashing and beeping signs of all process steps.
- To avoid ambience heating and save energy, an alarm will sound after 5 minutes if the equipment is not being operated and a "Prot" message will also be displayed warning
- that the muffle should be partially and manually closed.
 - Night Mode: after 2 hours without any operation, the muffle is closed and the temperature is kept at 100°C.
 -

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
- Operation time totalizer showed in hours.
- Low cost muffle replacing system.
- Elevator with bottom entrance perfectly balanced
- Electronic system forced ventilation keeps the components safe.
- ARGÔN SISTEM system, which allows perfect cleaning and helps with the ceramic burn accentuating the colors, enhancing both the translucence and opacity.
- Serial communication port for check-up (technical assistance)
- Safety temperature devices, thermocouple, vacuum, elevator, power, piston and compressed air system.
- Liquid crystal panel with 4 lines by 20 characters on a blue background.
- Touch-sensitive thin layer with sound alarm recognition for the operator's touch.

3 Installing

The numbers in parentheses refer to the pictures that follow.

Your Alumni-Press-2 should be far from curtains and inflammable material. A furnace is a heat generator which needs to be dissipated; otherwise its components will be overheated. Therefore, place the equipment at a ventilated area which allows free air circulation. It is advisable to leave a 15cm clearance between the furnace and anything that can impair the ventilation. Place the furnace far from taps or basins which can sprinkle water over the equipment.

PIC.1

- (1) - Network voltage identifier / power cable
- (2)- Fuse holder
- (3)- Main switch on / off
- (4)- Control system forced ventilation
- (5)- Socket for plugging the vacuum hose
- (6)- Socket for plugging the vacuum pump
- (7)- Argon system hose connection
- (8)- Serial communication
- (9)- Muffle ventilation
- (10)- Compressed air hose connection
- (11)- Press system pressure regulator
- (12)- Press system pressure manometer

Notes

- 3-1) Check if your network voltage is the same as the one shown on the furnace label Pic 1 – item (1).
- 3-2) Install the furnace into an exclusive electrical network using 6mm wiring to a 110 volts network, or 4mm wiring to a 220 volts network.
- 3-3) Never connect the furnace to the same network in which other furnaces, compressors, electrical taps or any other high electricity consumption devices are already connected.
- 3-4) Check if the outlet in which the furnace will be plugged in is in perfect condition and of good quality to bear at least 20 amperes.
- 3-5) Connect the grounding end (round pin) to a grounding rod and never to the electrical network neutral.
NOTE: Not checking the procedures above will interfere on the equipment correct working, on its guarantee and also on the operator safety.
- 3-6) If your electrical network is very instable (the lights flash a lot), it is advisable to use a voltage regulator such as a saturated autotransformer with at least 2kw capacity. Never use, by any means, regulators used in computers.
- 3-7) Connect the vacuum hose to the furnace Pic-1 item (5) Plug the vacuum pump power cable to the outlet which is at the backside of the furnace Pic-1 item (6)
- 3-8) Connect the compressed air hose to the furnace. Pic-1 item (10). The air should be supplied dried by the compressor at minimum pressure of 2 kg/cm² (30lbf/pol²). This pressure should be kept while the furnace is working with pressed ceramic. The lack of pressure during this process will cause the process interruption and a message will be shown on the control display.
- 3-9) Connect the argon hose to the furnace and to the cylinder outflow regulator Pic-1 item (7). The argon used is the commercial one. The cylinder and the outflow regulator are not supplied with the equipment, we recommend you to contact a gas supplier in your city. The maximum argon outflow should be adjusted at 10 L/m (liters per minute). Note.: We remind you that when using argon, both translucence and opacity are enhanced, but its use can be dispensed. To do so, just disconnect the correspondent hose and do not program its function. Check the hoses connections if they are very tight to avoid undesirable leaks.


15-18) KEYS (.) PONTO, 0,1,2,3,4,5,6,7,8,9 :- These keys allow the insertion of data for the recipes variables, activate the desired key and insert the new value.

15-19) KEYS WITH SYMBOLS:- These keys allow the insertion of personalized identification for the recipe names at the data editing and saving moments.

Elevator times table.

ELEVADOR				TEMPO APROXIMADO	
S	M	.	D	M	Sai do fim de curso
S	D	.	D	D	Percurso contínuo
S	1	.	D	1	1 Minuto
S	2	.	D	2	2 Minutos
S	3	.	D	3	3 Minutos
S	4	.	D	4	4 Minutos
S	5	.	D	5	5 Minutos
S	6	.	D	6	6 Minutos
S	7	.	D	7	7 Minutos
S	8	.	D	8	8 Minutos
S	9	.	D	9	9 Minutos

KEYS

16 Specifications


- Control - Micro processed state-of-the-art MOTOROLA.
- SMD electronic circuit cards
- Liquid crystal display with BACKLIGHT. Blue background.
- Power = 1800 Watts.
- Power supply 110v. or 220v. 50/60 Hz. Under request.
- Fuses: 20 amperes to 110 volts, 15 amperes to 220volts.
- Communication RS-232 – for technical assistance.
- Alphanumeric keyboard with 39 programming and monitoring keys.
- Maximum operation temperature 1.180°C
- Nominal working pressure 60 Lbf / in2 (Piston)
- Nominal working vacuum –27 in Hg. / mm Hg.
- Furnace dimensions:
- Width = 330 mm
- height = 800 mm
- depth = 420 mm
- net weight = 22,5 Kg.
- Gross weight = 28,5 Kg.
- Packaging dimensions 440mm x 840mm x 490 mm

4 Control panel

4-1) The control panel has a liquid crystal display with 4 lines by 20 columns on a blue background with yellow printing. It also has inner BACKLIGHT illumination, which allows easy identification of the programmed functions and register monitoring.

4-2) It has a touch-sensitive thin layer which facilitates the operator's job when activating its functions which are also signaled by beeping sounds.

4-3) Panel identification and its main functions. PIC-2

4-4) DISPLAY - Liquid crystal display 4 / 20 blue. Item-21.

4-5) PROGRAM – Sequential program screen. Item-22.

4-6) MONITORING – Control system monitoring key. Item-23.

4-7) QUICK FUNCTIONS – Use the PROGRAM key, and sequential screens for programming and monitoring. Item-24.

4-8) ELEVATOR ASCENT AND DESCENT – digital command turns the elevator motor on and off. Item-25.

4-9) AUXILIARY FUNCTIONS – Argon washing, quick cooling, pre burn, pressing. Item-26.

4-10) NIGHT – Automatic or manual resting system, used only when operating with ceramic. Item-27.

4-11) AUTO-START – This function allows the burn start in automatic position starting from ambient temperature. Item-28.

4-12) START/STOP - This function allows you to start or stop the burn process through command. Item-29.

4-13) KEYS SETS WITH GRAY BACKGROUND COLOR – Exclusive operation of the control and programming systems. Item-30.

4-14) NUMERIC KEYBOARD –Keyboard used to add numeric content to the recipes variables. Item-31.

4-15) MESSAGES – Field destined to the process stages visualization. Item-32.

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5 System start up - step by step

5-1) Plug the furnace into an exclusive and of good quality socket which should be correctly grounded to a copper rod and never to the electrical network neutral. Check if the network voltage is the same as the one in the equipment, which is shown on the power cable label or printed at the furnace backside. Pic-1 item (01).

5-2) Turn the main switch on Pic-1 item (3) and wait for the equipment to start-up. The first screen to be shown has the data for the programs versions. Ex-1.

STATUS	REC°N	TEMP °C	VÁCUO
E D G	E Q U I P A M E N T O S		
	E D G H A N D V 2 . 2		
	E D G B A S E V 1 . 3		
W W W . E D G . C O M . B R			
ELEVADOR	MENSAGENS		

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5-3) The following screen shows the initializations and check-up for the entire system and a scroll will keep running across the screen while you perform the initial configurations. Ex-2.

STATUS	REC°N	TEMP °C	VÁCUO
I N I C I A L I Z A N D O			
		◀ ▶	
ELEVADOR	MENSAGENS		

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5-4) After the checking-up is done, the system loads the ceramic standard no. 01 recipe and will keep waiting for the next instruction by the operator. We remind you this recipe will be the base for new recipes and the program does not allow it to be erased. Ex-30

STATUS	REC°N	TEMP °C	VÁCUO
M O N I T	0 0 1	0 0 2 5	0 0
T 1 =	0 6 0 0		
V E L 1 =	0 6 0		
S 1 D D . P A D R ã O C E R A M . .			
ELEVADOR	MENSAGENS		

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5-5) The furnace must not start the heating process without the heating spool put on the tray, otherwise the sealing ring and the tray will be seriously damaged. Use adequate accessories according to the wanted procedures PIC-4 e PIC-5.

15-7) KEY - GRAVA:- Activate the GRAVA key to save into the control memory a new recipe. Choose the number, the name and activate the GRAVA key once again to confirm the operation.

15-8) SETAS :- Use the right and left arrows to scroll among the selected digits during the programming stage and also to activate the argon fields filling, cooling and auto start.

15-9) KEY - OK:- This key is used to scroll among the selected digits during the programming stage and also to confirm the receipting of messages in general.

15-10) KEY - ENVIA : This key is used when we want to send a new recipe to the furnace manually or when we alter just one datum and we want it to be performed. We remind you that the ENVIA function occurs automatically every time that we move to the monitoring screen.

15-11) KEY - START STOP:- This key has double function, when it is activated with the muffle opened it will be interpreted as START, when activated with the muffle closed it will be interpreted as STOP.

15-12) KEY - ARGON:- This function allows the washing with argon automatically at the beginning of the burn process, activate the PROGRAMA key until you reach (tela 4) and confirm the ARGON function with the LEFT ARROW filling the relative field.

15-13) KEY - RESFR:- This function helps with the forced cooling between burn processes, using the vacuum pump to circulate cold air in the muffle helping with the heat exchange. Activate the PROGRAMA key until it reaches (tela 4) and confirm the ARGON function with the LEFT ARROW filling the relative field.

15-14) KEY - NIGHT:- This key allows that the furnace is kept in a resting state while it is not being used. Ex. at night, it keeps the muffle semi closed at a temperature of 100°C.

15-15) KEY - AUTO START:- It allows the beginning of a new burn process without having to wait until the temperature reaches T1. Ex. After turning the furnace on, program it, place the piece to be burned and activate AUTO START. The process will be performed automatically.

15-16) KEY - SOBE:- It allows the electric elevator activation at the SOBE function with the operator touching. When activated it will remain on until the operator touches it again to interrupt the operation. It is turned off automatically when it reaches the end of high course. To program drying values, use the PROGRAMA key (tela-5) and type the desired value, check the times table below.

15-17) KEY - DESCE:- It allows the electric elevator activation at the SOBE function with the operator touching. When activated it will remain on until the operator touches it again to interrupt the operation. It is turned off automatically when it reaches the end of low course. To program slow exit values, use the PROGRAMA key (tela-5) and type the desired value, check the times table below.


5-8) Move the elevator movable tray system to the front item-(13) using the front knob and place the spool item-(14) on its mortise, next take the tray back to its correct position on the base to start the equipment.

6 Programming conventional ceramic recipes

6-1) For better explanation, we will follow practical examples for conventional and pressed ceramic recipes.

EX-1 RECIPE EXAMPLE - CONVENTIONAL CERAMIC

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6-2) Start the programming reading the Standard recipe that comes from the factory (01), it will serve as reference to program the further recipes.

6-3) Activate the LER key, choose the recipe 01 typing number 1 in the REC field, or scroll through the options list using the left or right arrows. The recipe name will be shown on the message field (PADRÃO CERAM) EX-5, at this point you should confirm your option activating the LER key again to complete the operation. The recipe will be automatically sent to the base and, when there is no recipes saved the message (RECEITA VAZIA) will be shown.

Errorr 2 – Elevator flaw. ** ELEV-SOBE **

A flaw on the elevator ascent procedure is interrupted indicates electromechanical problems. Check if there are no obstacles blocking the elevator movement. Check - OK

Errorr 3 - Elevator flaw ** ELEV-DESCE **

A flaw on the elevator descent procedure is interrupted indicates electromechanical problems. Proceed as the previous item. Check - OK

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Indicates that the speed is adjusted at zero at T2. Check - OK

Errorr 5- “heating speed equals zero” ** VELOCIDADE **

Indicates that the speed is adjusted at zero at T1. Check - OK

Errorr 6 –Attempt to change the standard recipe** REC-PADRÃO **

The standard recipes cannot be changed. Save the new recipe in another available position in the memory. Check - OK

ERROR - 07 - Vacuum flaw – ** FALHA-VÁCUO **

Message:- It was not possible to crate the vacuum. Check - OK

ERROR - 08 - Vacuum flaw – ** FALHA-VÁCUO **

Message:- It was not possible to turn the vacuum off. Check - OK

ERROR - 09 - Version updating flaw ** PROGRAMA **

Message:- Current furnace base incompatible version. Check - OK

ERROR - 10 - A/D conversion flaw ** TEMPERATURA **

Message:- Temperature over the allowed limit. Check - OK

ERROR - 11 - Piston flaw ** COMPRESSOR **

Message:- The compressed air pressure is too low. Check - OK

ERROR - 12 - Programming flaw ** PRENSAGEM **

Message:- The pressing time should be shorter that the burn one. Check - OK

ERROR - 13 - Programming flaw ** PRENSAGEM **

Message:- Vacuum and burn time should be teh same Check - OK

ERROR - 14 - Programming flaw ** ARGÔNIO **

Message:- It is not possible to program argon washing without vacuum. Check - OK

ERROR - 15 - Programming flaw ** NIGHT **

Message:- This function can only occur in baseline T1 with the elevator on the base. Check – OK

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Message:- It was not possible to connect the furnace. Check - OK

ERROR - 17 - Movable tray device flaw ** START **

Message:- Elevator tray out of the correct base. Check - OK

ERROR - 18 - Technical access flaw ** ROT-TÉCNICA **

Message:- Access only allowed to authorized personnel. Check - OK


ERROR - 19 - Operator access flaw ** SENHA **

Message:- The password does not match, type the password again. Check – OK

15 The keyboard and its functions

15-1)The alphanumeric keyboard has 39 programming functions and burn system monitoring keys, divided in blocks, numeric, alphabetic, auxiliary functions and control.

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
15-3) KEY-MONITORAÇÃO:- This key allows the burn variables visualization during the process. On a single screen the data are updated according to the stages being performed.


15-4) KEY - A/N:- This key allows altering the writing mode between numeric and alphabetic according to the programming necessities when editing recipe names.

15-5) KEY – LIMPA REC :- This key can be used when editing and programming recipes, when we want to delete all the content of a screen, EX- T2, PAT2, VAC2, VEL2. or when we want to delete permanently a recipe from the memory. This function is only available in saving mode.


15-6) KEY – LER:- Activate the LER key to load a previously saved in the control memory recipe, type the LER key and the operation will be cancelled.

PIC.4

PIC.5

PIC.3

PIC. 4

5-6) Accessories for CONVENTIONAL CERAMIC burn.

A) Use the spool – item (14)

B) Use the rigid mat on the spool. item (15)


C) Use isothermal pins for supporting the piece Item (16).

PIC.5

5-7) Accessories for PRESSED CERAMIC burn.

D) Use the spool – item (17).

E) Use the refractory platform item (18) in the spool to accommodate the ring.

PIC.6

Warning: The set of RINGS, TABLETS AND PLUNGERS, should not surpass the maximum height of 92 millimeters, as shown on PIC-(6).

13 Alarms and messages.

13-1) To avoid damage to the equipment and programming mistakes, several safety devices work in the flaws prevention system. Messages and Error codes will be visualized on the display.

13-2) Operation messages :- All messages being performed are shown on the field reserved for the messages on the screen. Each message will be shown at its time according to the current burn stage, for example:- RAMPA, PATAMAR, PRÉ QUEIMA, RESFRIAMENTO, AUTO START, among other, these messages are intended to illustrate the burn process being performed.

13-3) Programming messages:- They accuse programming flaws which are being saved or sent to the furnace, and they do not allow that wrong information, which can damage the burn, is performed.

13-4) System flaw messages:- Electrical problems may occur and can interfere with the equipment perfect working. To avoid the flaws, some vital functions are monitored continually and are communicated if there is some problem.

14 Messages definitions

PROTECTION MODE :- It protects the equipment against high temperatures on the panel and also avoids energy waste. It is activated after 5 minutes with no operation after the elevator is back to the base.

14-2) NIGHT MODE: - This function puts the furnace in resting state and also at low temperature closing the muffle, it can be used at night.

14-3) ERRORS MODE: - These are messages that communicate general problems and they will stay on screen until the operator recognizes them and take the correct procedures and confirm them using the OK key.

Error 1- "Thermocouple flaw" ** TERMOPAR **.

The process being performed is aborted, an alarm will sound warning about the flaw. Contact the technical assistance. Check - OK

STATUS	REC°N	TEMP °C	VÁCUO
L E R			
R E C : 0 0 1		C E R Â M I C A	
A L F A			
N O M E : . P A D R Ã O . C E R Â M .			
ELEVADOR	MENSAGENS		

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6-4) After the standard ceramic recipe is read, the system returns to the monitoring screen automatically. EX- 11

STATUS	REC°N	TEMP °C	VÁCUO
M O N I T	0 0 1	0 6 0 0	0 0
T 1 = 0 6 0 0			
V E L 1 = 0 6 0			
S D . D D P A D R Ã O . C E R Â M			
ELEVADOR	MENSAGENS		

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6-5) To start programming a new recipe, activate the PROGRAMA key and the display will show the message TELA1, the cursor will be flashing over the field that refers to the temperature T1=600. At this stage, we have two options. We can simply type the new values on the fields which have been already filled or clean up all the values by using the LIMPA REC key for typing in the empty fields. EX-12.

Note: For CONVENTIONAL CERAMIC recipe EX-1, we will use the empty fields activating the LIMPA REC function.

STATUS	REC°N	TEMP °C	VÁCUO
T E L A 1	0 0 1		
T 1 = 0 0 0 0			
N O M E :			
ELEVADOR	MENSAGENS		

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6-6) On TELA-1, with the fields empty, the cursor will be flashing over the field T1=0000. Type the new value T1=0600 following the conventional ceramic model EX-6 step-by-step.

STATUS	REC°N	TEMP °C	VÁCUO
T E L A 1	0 0 1		
T 1 = 0 6 0 0			
N O M E : P A D R Ã O . C E R Â M			
ELEVADOR	MENSAGENS		

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Activate the PROGRAMA key again and a new screen TELA-2 will be shown, and the cursor will be flashing over the field T2=000. Type the new value T2=960 and confirm using the OK key. Next, the cursor will flash over the velocity field VEL2=000, type the new value VEL2=55 and confirm using the OK key.

Now type the new value for baseline time PAT2=01:00 and vacuum time VAC2=01:00 these times are programmed in five minutes. Confirm the operation using the OK key EX-7.

STATUS	REC°N	TEMP °C	VÁCUO
TELA 2	0 0 1		
T 2 =	0 9 6 0	P A T 2 =	0 1 : 0 0
V E L 2 =	0 5 5	V A C 2 =	0 1 : 0 0
N O M E :	P A D R ã O . C E R ã M		
ELEVADOR	MENSAGENS		

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Activate the PROGRAMA key again and a new screen TELA-3 will be shown, and the cursor will be flashing over the field T3=000. Type the new value T3=980 and confirm using the OK key. Next, the cursor will flash over the velocity field VEL3=000, type the new value VEL3=020 and confirm using the OK key. Now type the new value for baseline time PAT3=01:00 and vacuum time VAC3=00:00 these times are programmed in five minutes. Confirm the operation using the OK key EX-8.

STATUS	REC°N	TEMP °C	VÁCUO
TELA 3	0 0 1		
T 3 =	0 9 8 0	P A T 3 =	0 1 : 0 0
V E L 3 =	0 2 0	V A C 3 =	0 0 : 0 0
N O M E :	P A D R ã O . C E R ã M		
ELEVADOR	MENSAGENS		

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6-9) Activate the PROGRAMA key again and a new screen TELA-4 will be shown, and the cursor will be flashing over the pre burn field PRÉQ=00:00, type the new value PRÉQ=05:00 and confirm the operation using the OK key. Next, the cursor will be flashing over the pressing field PREN=00:00, keep the value at ZERO and confirm the operation using the OK key. Now the cursor will flash over the argon field ARGON, to confirm this option activate the left arrow and the correspondent field will be filled, confirm the operation using the OK key, ARGON=SIM ATIVADO, repeat the previous step to program the quick cooling field RESFR=NÃO and also for the auto start field AUTOS=NÃO. When these fields are filled, the selected function will be activated and the parameters will be part of the recipe. EX-9.

STATUS	REC°N	TEMP °C	VÁCUO
TELA 4	0 0 1		
P R É Q =	0 5 : 0 0		
P R E N =	0 0 : 0 0		
N O M E :	P A D R ã O . C E R ã M		
ELEVADOR	MENSAGENS		

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12-2) ARGON SISTEM:- This function performs the automatic washing of the piece at the beginning of the burn process. Use only inert and not oxidizing argon gas, it will help contrasting the colors, enhancing the translucence and opacity. EX-55.

12-3) QUICK COOLING:- This function is used to optimize the time spent between the burn processes, forcing the heat exchange so that the temperature can return as quick as possible to T1 and allowing the beginning of a new burn process EX-55.

12-4) NIGHT:- This functions is automatically activated every time the equipment is idle for two hours time, or it can be activated manually using the NIGHT key, which will put the furnace under total resting. When that happens the elevator is retracted closing the muffle, then the temperature decreases from T1 to 100 ° C. This option can be activated only with baseline T1 and the elevator on the base for ceramic burn. The main objective of this function is to save energy, avoiding the waste and keeping the muffle heated at a minimum temperature of 100°C, avoiding moisture. NOTE: - To deactivate this function, activate the NIGHT key once again.

12-5) AUTO START:- This option allows the automatic burn start when the temperature reaches the T1 baseline with the elevator on the base, at this moment the burn will be performed. To do so, it is necessary that the piece is already placed on the burn tray. A visual and sound alarm will start communicating the process beginning. NOTE: The burn will be performed once.

This option can be activated both from the T1 ascent ramp from 25° C to 600 ° C and descent from 960° C to 600° C at cooling, directly from the AUTO-START key, or even be programmed using the PROGRAMA key - (tela 4), in which it will become part of the recipe. EX-55

Programming auxiliary functions.

STATUS	REC°N	TEMP °C	VÁCUO
TELA 4	0 6 5		
P R É Q =	0 5 : 0 0		
P R E N =	0 0 : 0 0		
N O M E :		
ELEVADOR	MENSAGENS		

PROGRAMANDO
PRÉ-QUEIMA
ARGON
RESFRIAMENTO
E AUTO-START
EX - 55

10-3). System flow monitoring. The equipment has several safety devices which allow good performances, if there is any kind of defect the burn will be automatically aborted and an Error message will be on the display until the operator fixes the problem and releases the operation using the OK key. Following we will give you some common Errors examples

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TELA DE MENSAGENS
ERRO DE START
BANDEJA-ELEV.
EX - 93

11) QUICK FUNCTIONS. These functions allow the programming or monitoring of the desired variables using the shortcut from the PROGRAMA key.

12) AUXILIARY FUNCTIONS

12-1) AUXILIARY FUNCTIONS:- They are complementary resources which can help during the burn processes improving the product final quality. It facilitates the control system operation protecting the equipment and it also minimizes the energy consumption speeding the process final time up. They are:- ARGON SISTEM, RESFRIAMENTO RÁPIDO, AUTO-START, PRÉ-QUEIMA, NIGHT, EX-55.

6-10) Activate the PROGRAMA key again, and a new key TELA-5 will be shown in which the cursor will be flashing over the “D” field for the SOBE function. Take the cursor to the position (5) using the left and right keys and confirm the operation using the OK key. The confirmation of this option will be shown in parentheses on the SOBE field. Repeat the previous step to program the elevator function DESCE, move the cursor to the position (3) and e confirm the operation using the OK key. EX-10.

STATUS	REC°N	TEMP °C	VÁCUO
TELA 5 0 0 1 ELEVADOR			
SOBE (5) MD 1 2 3 4 5 6 7 8 9			
DESCE (3) MD 1 2 3 4 5 6 7 8 9			
NOME : . . . PRENSADA . . .			
ELEVADOR	MENSAGENS		

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
6-11) Activate the PROGRAMA key again, the system will return to TELA-1 showing the programming beginning. At this point the recipe is ready to be saved GRAVADA in the memory, check the step-by-step for saving recipes PG-18-item-9. To send the recipe to the base just activate the MONITORAÇÃO screen.

7 Programming pressed ceramic recipes

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7-1) Start the programming reading the Standard recipe that comes from the factory (01), it will serve as reference to program the further recipes.

7-2) Activate the LER key, choose the recipe 01 typing number 1 in the REC field, or scroll through the options list using the left or right arrows. The recipe name will be shown on the message field (PADRÃO CERAM) EX-5, at this point you should confirm your option activating the LER key again to complete the operation. The recipe will be automatically sent to the base and, when there is no recipes saved the message (RECEITA VAZIA) will be shown.

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7-3) After the standard ceramic recipe is read, the system returns to the monitoring screen automatically. EX- 11

STATUS	REC°N	TEMP °C	VÁCUO
MONIT	0 0 1	0 6 0 0	0 0
T 1 = 0 6 0 0			
VEL 1 = 0 6 0			
S D . D D . PADRÃO . CERÂM			
ELEVADOR	MENSAGENS		

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7-4) To start programming a new recipe, activate the PROGRAMA key and the display will show the message TELA1, the cursor will be flashing over the field that refers to the temperature T1=600. At this stage, we have two options. We can simply type the new values on the fields which have been already filled or clean up all the values by using the LIMPA REC key for typing in the empty fields. EX-12.

STATUS	REC°N	TEMP °C	VÁCUO
TELA 1	0 0 1		
T 1 = 0 0 0 0			
NOME :			
ELEVADOR	MENSAGENS		

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NOTE: For CONVENTIONAL CERAMIC recipe EX-2, we will use the empty fields activating the LIMPA REC function.

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9-3) After the ceramic recipe is saved, the system returns to the monitoring screen automatically. EX- 30.

STATUS	REC°N	TEMP °C	VÁCUO
MONIT	0 0 2	0 5 8 0	0 0
T 1 = 0 6 0 0			
VEL 1 = 0 6 0			
S D . D D OPAÇO			
ELEVADOR	MENSAGENS		

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10 System monitoring

10-1) Burn processes monitoring. This screen allows the visualization of all variables for the burn during the process. Each burn stage will be shown on the message field of the display, making it easier the operation monitoring. Process messages

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MENSAGENS			
L	A	V	A
R	E	S	F
I	N	I	C
P	R	E	N
A	U	T	O
R	E	S	F
S	T	O	P
D	E	S	C
P	A	D	R
S	U	B	S
N	O	M	E

10-2) The variables like T1, T2, T3, Vel1, Vel2, Vel3, Pat1, Pat2, Pat3, Vac2, Vac3, Pre-burn, Pressing, recipe number, elevator sobe (D) and desce(D), will be continually updated during the burn process on the monitoring screen. EX-38.

NOTE: The saved recipes names will be shown on the first burn, ambient temperature burn until it reaches the baseline T1.

Items in contrast are part of this screen EX-38

STATUS	REC°N	TEMP °C	VÁCUO
MONIT	0 0 2	0 6 0 0	- 2 7
T 2 = 0 9 6 0		PAT 2 = 0 3 : 0 0	
VEL 2 = 0 4 0		VAC 2 = 0 2 : 0 0	
S D . D D RAMP A			
ELEVADOR	MENSAGENS		

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STATUS	REC°N	TEMP °C	VÁCUO
L E R			
R E C : 0 0 3		C E R Â M I C A	
N O M E : . . .		P R E N S A D A	
ELEVADOR	MENSAGENS		

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8-5) After the standard ceramic recipe is read, the system returns to the monitoring screen automatically. EX- 31
Pressed ceramic monitoring REC-03

STATUS	REC°N	TEMP °C	VÁCUO
M O N I T	0 0 3	0 8 0 0	0 0
T 1 = 0 8 0 0			
V E L 1 = 0 6 0			
S D D D . . .		P R E N S A D A	
ELEVADOR	MENSAGENS		

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9 SAVING RECIPES

9-1) SAVING RECIPES IN THE DATA BASE.

Ex-35 – Conventional ceramic opaque burn. REC-02

9-2) EX-35 - Opaque burn. Activate the GRAVA key, choose recipe 02 typing number 02 on the REC field, or scroll through the options list using the left or right arrows. The recipe name will be shown on the message field [. OPACO.]. If you wish to overwrite the pré saved recipe, confirm the operation activating the GRAVA key, in this case a message will be shown [SUBST. RECEITA? (S/N)] if you want to choose yes SIM type the letter [S] the operation will be concluded and the recipe will be sent to the furnace base. If you want to choose no NÃO type the letter [N] the cursor will return to the REC=002 field so that the recipe number can be altered. To change the recipe name, activate the OK key, the cursor will move to the field referring to the recipe, type a new name, EX. (. . . OPACO. . .) confirm the option activating the GRAVA key again. The recipe will be automatically sent to the furnace base. EX-35.

Saving opaque burn. REC-02

7-5) On TELA-1, with the fields empty, the cursor will be flashing over the field T1=0000. Type the new value T1=0800 following the pressed ceramic model EX-13 steps.

STATUS	REC°N	TEMP °C	VÁCUO
T E L A 1	0 0 1		
T 1 = 0 8 0 0			
N O M E : . . .		P R E N S A D A	
ELEVADOR	MENSAGENS		

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7-6) Activate the PROGRAMA key again and a new screen TELA-2 will be shown, and the cursor will be flashing over the field T2=000. Type the new value T2=960 and confirm using the OK key. Next, the cursor will flash over the velocity field VEL2=000, type the new value VEL2=55 and confirm using the OK key. Now type the new value for baseline time PAT2=01:00 and vacuum time VAC2=01:00 these times are programmed in five minutes. Confirm the operation using the OK key EX-14.

STATUS	REC°N	TEMP °C	VÁCUO
T E L A 2	0 0 1		
T 2 = 0 9 6 0		P A T 2 = 2 9 : 0 0	
V E L 2 = 0 5 5		V A C 2 = 2 9 : 0 0	
N O M E : . . .		P R E N S A D A	
ELEVADOR	MENSAGENS		

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7-7) Activate the PROGRAMA key again and a new screen TELA-3 will be shown, and the cursor will be flashing over the field T3=000. For this recipe the stage three (três) will not be used, so all the fields for this screen should remain ZERO. EX-15
T3=000, VEL3=000, PAT3=00:00, VAC3=00:00

STATUS	REC°N	TEMP °C	VÁCUO
T E L A 3	0 0 1		
T 3 = 0 0 0 0		P A T 3 = 0 0 : 0 0	
V E L 3 = 0 0 0		V A C 3 = 0 0 : 0 0	
N O M E : . . .		P R E N S A D A	
ELEVADOR	MENSAGENS		

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7-8) Activate the PROGRAMA key again and a new screen TELA-4 will be shown, and the cursor will be flashing over the pre burn field PRÉQ=00:00, keep the values at ZERO and confirm it using the OK key. Next, the cursor will be flashing over the pressing field PREN=00:00, type the new value PREN=09:00 and confirm the operation using the OK key. Now the cursor will flash over the argon field ARGON, keep it deactivated ARGON=NÃO and confirm the operation using the OK key. To program the quick cooling field, activate the left arrow to fill the field RESF=SIM, confirm the operation using the OK key. The cursor will flash over the autostart field, keep it deactivated and confirm the operation using the OK key AUTOS=NÃO. When these fields are filled, the selected function will be activated and the parameters will be part of the recipe. EX-16

STATUS	REC°N	TEMP °C	VÁCUO
TELA 4	0 0 1		ARGON
PRÉQ = 00 : 00			RESFR
PREN = 09 : 00			AUTOS
NOME : . . .	PRENSADA		
ELEVADOR	MENSAGENS		

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7-9) Activate the PROGRAMA key again, and a new key TELA-5 will be shown in which the cursor will be flashing over the "D" field for the SOBE function. Take the cursor to the position (5) using the left and right keys and confirm the operation using the OK key. The confirmation of this option will be shown in parentheses on the SOBE field. Repeat the previous step to program the elevator function DESCE, move the cursor to the position (D) and e confirm the operation using the OK key. EX-17.

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7-10) Activate the PROGRAMA key again, the system will return to TELA-1 showing the programming beginning. At this point the recipe is ready to be saved GRAVADA in the memory, check the step-by-step for saving recipes PG-18-item-9. To send the recipe to the base just activate the MONITORAÇÃO screen.

8 Reading recipes

8-1) READING RECIPES FROM THE DATA BASE.

Ex-27 – Opaque burn for conventional ceramic. REC-02

Ex-28 – Burn for pressed ceramic. REC-03

8-2) EX-27 - Opaque burn. Activate the LER key, choose the recipe 02 typing number 2 in the REC field, or scroll through the options list using the left or right arrows. The recipe name will be shown on the message field (. . . OPACO. . .), at this point you should confirm your option activating the LER key again to complete the operation. The recipe will be automatically sent to the furnace base EX-27. Opaque burn, conventional ceramic. REC-02

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8-3) After the standard ceramic recipe is read, the system returns to the monitoring screen automatically. EX- 30

STATUS	REC°N	TEMP °C	VÁCUO
MONIT	0 0 2	0 5 8 0	0 0
T 1 = 0 6 0 0			
VEL 1 = 0 6 0			
S D D D . . .	OPACO		
ELEVADOR	MENSAGENS		

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8-4) EX-28- Pressed ceramic burn. Activate the LER key, choose the recipe 03 typing number 3 in the REC field, or scroll through the options list using the left or right arrows. The recipe name will be shown on the message field (. . . PRENSADA. . .), at this point you should confirm your option activating the LER key again to complete the operation. The recipe will be automatically sent to the furnace base EX-29. Pressed ceramic burn. REC-03