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powerCAST *1700*

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		20
--	--	----

Número de Série

--

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--

1 Safety recommendations

Warning: Read this manual carefully before installing and operating the equipment. The powercast is a high frequency centrifuge, projected to cast and centrifuge all the alloys used in dentistry and jewelry store, except for titanium. Any use that is not in accordance with the one that was projected to will be considered as improper use, as well as non compliance of the rules of operation, safety and maintenance indicated in this manual.

2 Receipt

When receiving the equipment, before signing the term of receipt, check the package general state; if there is any damage contact the factory immediately and you will have procedure orientation. We remind you that the good travels by the buyer's own responsibility and it is secured by the transportation service.

After opening the package check its contents:

- Equipment.
- 3 ceramic melting pots.
- 3 ceramic receptors.

- 1 water entrance hose.
 - 1 water exit hose.
 - 1 steel receptacle.
 - 1 water pressure inspection device.
 - 1 bracket.
 - Instruction manual.
 - Warranty certificate.
 - Equipment with argon device is accompanied by one more hose and 2 brackets.
- Description of the powercast induction centrifuge

3 Operation

The powercast is a high frequency induction centrifuge, strong and easy to operate, proper to cast and centrifuge all the dental alloys and jewelry store, except for titanium and its alloys.

Its compact design occupies the same space of a conventional centrifuge. It produces high frequency magnetic fields that induced in the metal cause its heating and consequent fusion. The magnetic field also promotes a strong agitation in the casted metal, homogenizing the alloy and improving its metallurgic features.

The heating by induction is clean, minimizing the formation of oxides and eliminating the carbon absorption, commonly verified with the use of soldering pipes.

The non-use of dangerous gases like LPG and oxygen provide to the powercast extreme safety and reliability giving a total casting control, warranting the repeatability of the results independent on the operator.

All the casting and centrifuge occur inside the closed chamber eliminating the possibility of accidents.

Your powercast can be equipped with the exclusive argon system® that consists of argon gas injection in the melting pot during the cast and centrifuge.

Argon is an inert gas that when injected to the metal expulses the oxygen, minimizing the oxidation during the casting. The process is totally automatic; the gas is injected at the beginning of the casting and it remains this way during the centrifuge. We also have rotation, acceleration and balance control.

The powercast is:

- Safe.
- Economic.
- It produces high quality homogeneous casts.

- It practically eliminates the necessity of reworking.
- Easy handling.
- It does not occupy much space..
- It reduces the waste of material.

Warning: This equipment generates electromagnetic waves. It must not be operated by individuals who use pacemakers. Do not introduce combustible materials that can cause fire inside the chamber during the operation. Do not look at the casted metal without using the dark protection filter. This equipment is not proper for the casting of light metals like aluminum, magnesium, titanium, which can cause dangerous situations.

4 Installation


Installation local

Your powercast must be installed on a leveled, strong, wood or brick-built base, between 0,65 and 0,75m height, which allows the operator to observe the casting process through the view finder without any difficult.

The local must be cleaned and aired, it is important to provide a distance of about 20cm from any revetment in order to allow the free circulation of air needed to ventilation. It is not advisable that wax elimination furnaces be installed on the same benchwork.

Electrical installation

The powercast is produced at 220 Volts, it consumes the maximum of 3.000 Watts, 14 amp. Provide a supply net with 2,5mm² wire, install a 20 amp. thermal circuit breaker. Preferably, use the services from a professional.

- (01) Power cable
- (02) General switch
- (03) Argon entrance
- (04) Auxiliary outlet

According to the rules and for your own safety it is necessary the connection of the outlet central pin to a proper grounding; do not use the neutral of the electric network.

Warning: use an exclusive outlet; do not use branch-socket devices.

4.1 Refrigeration water

The induction coil and the high frequency generation system are refrigerated by circulating water, the good working of the equipment and its useful life depend on the observation of the following information:

The powercast needs a minimum water pressure of $1\text{kg}/\text{cm}^2$ ($15\text{lb}/\text{inch}^2$). The water must be clean and constant. If in your region the water has problems of sediments, it is necessary the installation of a proper filter. Otherwise, the system will be blocked causing damage to the equipment.

Use preferably water from the streets, which normally has higher pressure.

Provide the installation of a $3/4$ inch tap, garden type, about 1 meter from the equipment.

To check if the pressure in your network is proper use an inspection device that comes with the equipment (the following figure).

Provide: a 1-liter bowl (soda, alcohol bottle, etc.) a watch that counts the seconds or chronometer.

1. Connect the device of the figure to a tap.
2. Open the tap and measure how much time water takes to fill in the 1 liter bowl.
3. Repeat the operation twice and take notes on the time value.
4. Compare with the table below.

TIME (minutes)	PRESSURE (kg/cm^2)
3	0,5
02:30	1
02:00	2
01:30	3
01:00	4

If the time value obtained is less than 2:30 your equipment can be turned on normally, otherwise it will be necessary a pressurization pump that will be connected to the auxiliary outlet (4) .

Contact the factory and we will be able to guide you for the best solution.

Values above 4Kg/cm² are not advisable; a pressure reductor like the one used in buildings must be provided.


Remove the checking device from the tap and keep it. It can be useful later.

- Connect the hose that accompanies the equipment to the tap and the water entrance (5) to the back of the equipment. Open the tap and correct possible leakage.

- Connect the other hose to the water exit of the equipment (6), use a bracket to avoid leakage.

- The water used must be thrown away in a sink or drain. If this is not possible, provide a hose with a minimum of 25 mm internal diameter.

- Introduce the exit hose inside the biggest diameter and place it in a way that there is not reflux and that allows the water to flow by gravity until the closest drain.

Warning:

- Do not strangle the water exit hose.
- Do not reduce the length of the exit hose.
- Do not connect the exit hose to other one of same diameter.
- The actions above will cause poor operation of the equipment.

4.2 Argón

If your powercast is equipped with argon system®, it is necessary an industrial type argon gas cylinder and a pressure regulator coupled to it. The gas provider company of your city will guide you about the proper regulator model. A 3 m³ cylinder is enough for many castings since there is not any leakage in the connections.

Connect the hose to the pressure regulator and to the tip indicated as "Argon Entrance(3)" in the back of the equipment. Use brackets to warranty the sealing of the system. Rotate the regulator gauntlet totally to the left for later adjustment.

You can operate your equipment normally even if you haven't provided the installation of the argon cylinder.

4.2 Initial procedures


In order to best know your equipment and its operation, we advise the user to take test castings using as ring the steel receptacle (7) that comes with the equipment.


Warning - Danger:

The steel receptacle must not have any oil or water residual; the casted metal in contact with these components can cause serious damage to the equipment.

The numbers in brackets refer to the following figures.

- (16) Melting pot
- (17) Melting case
- (18) Ring support gauntlet
- (19) Balance gauntlet
- (20) Ring


- Open the door and remove the shims and other transportation accessories.
- Turn down the coil activating the releasing-cam (9) to the right.
- Put a melting pot (16) in its case (17).
- Put in the melting pot about 20 grams of metal, preferably nickel chrome. Distribute the tablets in the melting pot in a compact way.
- Rotating the ring support gauntlet (18) adjust the steel receptacle in a way that the orifice be aligned to the melting pot orifice.
- Tighten the ring support gauntlet until the receptacle is firm in its position.
- Loose the balance gauntlet (19) until the Centrifuge arm (15) works as a scale.
- Rotate the counterweight (22) until the arm is in the horizontal position.
- Tighten the gauntlet again, just touching it; it is not necessary to tighten it too much.
- Put the ceramic receptor (23) in its position.


Warning: Do not forget to use the ceramic receptor. If there is breakage of the melting pot the casted metal will be contained in the receptor not causing damage to the equipment. It is normal the receptor breaks when the metal contained in it is removed.

(22) Counterweig


(23) Ceramic receptor

(24) Argon injector

- If your powercast is equipped with argon system® put the argon injector (24) in the center of the melting pot.
- With the right hand rotate the Centrifuge arm (5) in a way that the bottom of the melting pot is on the coil (11), with the left turn up the coil lever (10).
- Adjust the centrifuge arm in a way that there is arm and coil lever locking; there is just one position where this occurs.

5 Control panel

- (25) Start key
Starts the casting process.
- (26) Acceleration switch
Adjust the start torque.
- (27) RPM switch
Adjust the number of rotations per minute of the centrifuge arm.
- (28) Preheating key
Preheating consists of heating the metal before putting the ring in its position maintaining its temperature (chapter 6-5).

(29) Temperature switch

Adjusts the shining of the reference filament according to the required temperature (chapter- 7) .

(30) Power switch

Adjusts the heating speedy. It can be activated during the fusion. It is very useful in the utilization of residuals when it is necessary to await all the metal pieces to reach the same temperature.

(31) Stop key

Stops the process at any time. It cancels the failure messages after the they are corrected.

As a guide, use the following tables:

Type of metal	Average density g/cm ³	Suggested Rotation rpm
Noble metals Ceramic gold	16 - 18	350
Yellow gold Low fusion	13 - 15	400
Seminoble	11 - 14	425
Nickel-chrome Chrome-cobalt Nonnoble metals	7 - 10	450

Metal weight gramas	Acceleration %
0 ~ 30	100
31 ~ 45	90
46 ~	80

6 Operation

Warning: For your own security, always use protection glasses.
Heated sparks can be thrown in the following operations.

- Turn on the General Switch (2) on the back of the equipment.
- Rotate the Temperature (29) switch totally to the left.
- Turn on the Service Switch (12).
- The display (32) will light up showing the version of the used software.

	W	W	W	.	E	D	G	.	C	O	M	.	B	R	
	I	N	D	U	C	A	O		V	1	.	0			

- And then the tests of water and refrigeration pressures will be taken.

	W	A	T	E	R		O	K							
	C	O	O	L	I	N	G		O	K					

- If there is any problem it will be shown as in the example:

			F	A	I	L	U	R	E	!	!				
					W	A	T	E	R						

Indication of lack of water or not enough pressure in the network.
Corrected the problem, press the STOP key and the process will be reset.

- If everything is ok the display will show the message:

					R	E	A	D	Y	!					
T	E	M	P		1	1	0	0	°	C					

If your powercast is not equipped with the temperature indicator, the message "TEMP" will not appear on the display. The temperature measure system will be detailed in the chapter 7.

- Place the RPM(27) switch in the position 450.
- Place the Acceleration (26) switch in the 100% position.
- Place the Power (30) switch in the "2" position.
- Close the lid (14) pressing it on the doorpost.
- If your powercast is equipped with the argon system® rotate the argon regulator gauntlet to the right until the exit indicator shows 10 liters/min. New adjustments

not necessary in the subsequent castings.

- Press the Start(25) key.
- The display(32) will show the message:

D	O	O	R											O	K
C	O	I	L											O	K

- If there is any problem the display will indicate it as the example:

				F	A	I	L	U	R	E	!	!			
				D	O	O	R		O	P	E	N			

Indicating that the door is open or closed improperly. Close the door and press the Stop (31) key and then Start key, the process will be reset.

Next the display will show the message below, indicating the beginning of the casting process.

M	E	L	T	I	N	G			0	:	0	0	M	I	N
T	E	M	P		1	1	0	0	°	C					

- The casting spent time will be shown on the upper right corner.
- Observing by the view finder (8) rotate the power switch to the right. The tablets will start being red and then they will be casted.
- Pull the dark protection view finder (32) at the moment that you find it necessary.
- When the metal is totally casted activate the releasing cam to the right.
- The centrifuge will be activated and the display will show the message:

				C	A	S	T	I	N	G					
T	E	M	P		1	1	1	0	°	C					

When the centrifugation is over wait for the total stop of the centrifuge;

- The display will show the message:

E	N	D			O	F		P	P	R	O	C	E	S	S
	M	E	L	T	I	N	G		0	:	2	5			

Indicating the casting time. This message will remain until the door opening.

Afterwards, the equipment will be ready for a new casting.

- Return the power switch for the "2" position.
- Open the door and using a clamp, remove the supporting receptacle and the casted metal of the receptacle.
- Cool the metal and the receptacle.
- With a spatula remove the residuals of the casting from the melting pot.
- Take several tests until you get totally familiarized with the process.
- In the real casting, the balance must be taken before putting the ring.

Warning: for your own security always use protection glasses. Heated sparks can be thrown in this operation. Remember that the receptacle must be perfectly dry before the use.

7 Real casting

For your own security keep the local work clean, use security glasses and proper tools to the work.

Rings

Any type of ring can be used, steel rings, free expansion system, PRP, not centered supply canal, etc.

You must be careful so that the canal of the ring metal entrance be aligned to the melting pot orifice.

It is necessary to attach the ring by means of the Ring Supporting Gauntlet.

It is not necessary to tighten the support on the ring too much.

Balance

When the system is unbalanced, strong turbulence is caused in the liquid metal causing erosion and bad filling in the modeling cavity.

That is why a good balance is important.

Make the balance with the rings before they go to the wax elimination furnace. Use a mark in the counterweight arm as reference.

Studies taken by us show that the heat loss in the center of the ring is about 10°C in the first minute, thus increase the temperature of the furnace in 10°C and make the thin balance with the ring in the casting position.

It is only necessary a small tighten in the balance gauntlet.

Procedure

Before removing the hot ring from the wax elimination furnace and begin the casting:

- Adjust the speed and acceleration parameters.
- Adjust the temperature. (equipment with measure system).
- Put the melting pot in its case.
- Put the metal to be casted in the melting pot.
- Put the centrifuge arm and raise the coil lever until its locking.
- Put the argon injector and check if the argon cylinder valve is open.

Casting

- Remove the furnace ring, put it correctly on its support and if you find it convenient use a proper ring mounting; even using the ring mounting is essential to tighten the ring supporting gauntlet. Take a thin balance.
- Close the door.
- Activate the start key.
- Rotate the power switch to the right.
- Use the view finder to observe the casting, wait for all the metal tablets reach the same temperature. Sometimes it is necessary to reduce the power so that there is not overheating of an only tablet.
- When you find it necessary use the dark glass view finder (35).
- When the metal is casted activate the releasing cam gauntlet.
- Wait for the total stop of the centrifuge.
- Open the door and remove the ring.
- The process can be interrupted at any time by pressing the "Stop" key or opening the door.
- At every new casting clean the melting pot. .

Warning: for your own security always use protection glasses; heated sparks from this operation can be thrown at you.

There is no limit for following castings. If, by chance, by any reason there is overheating the process will be interrupted and the display will show this message:

			F	A	I	L	U	R	E	!	!			
		O	V	E	R	H	E	A	T	I	N	G		

Press the "Stop" key, wait for 15 minutes and restart the process.

Preheating

To minimize the heating loss of the ring during the metal casting time the powercast has the preheating function that consists of preheating the metal in the melting pot before removing the ring from the furnace. To enable the function, activate the preheating key (28) on the panel before the beginning of the casting.

The display will show the message:

					R	E	A	D	Y	!					
T	E	M	P		1	1	1	0	°	C			P	R	E

Do not put the ring in its support, close the door, activate the Start key, adjust the power. The metal will be heated during 15 seconds, after this the heating will be interrupted and the display will show the message:

				W	A	I	T	I	N	G	!				
T	E	M	P		1	3	5	0	°	C			P	R	E


Open the door, remove the ring and put it in its support. Make the balance, close the door, the casting will follow regularly. It is not necessary to press the Start key. To disable this function press the preheating key again.

8 Temperature measure system

The operation of the temperature measure system is based on the principle that when a material is heated its temperature is proportional to the light emitted by it. In the powercast temperature measure system there is a filament that, by the passage of an electrical current gets incandescent to the temperature indicated on the display. When the casted metal has the same color of the filament both will be

in the same temperature.

To operate the system adjust the required temperature by rotating the "Temperature" (29) key on the panel, reading the value on the display.

When pressing the key "Start" the filament(33) lights up getting in the temperature selected on the panel.

Using the cursor(34) put the filament on the metal that is being casted.

When the filament is not visible and it starts to get dark the metal is in the temperature indicated on the panel. Activate the releasing cam.

When you find it necessary, use the dark protection view finder (35).

Take several testes to get familiarized with the system.

9 Melting pots

The melting pots are ceramics pieces with special shape and composition projected to support high thermal shocks that will be subjected to during the casting.

Its useful life depends on several factors among them: composition of the alloys that is used, casting temperature of the alloy, weight of the metal to be casted, etc.

Its dimensions can vary and sometimes it is necessary some settings so that it can have a perfect fitting in its case. Use a thick sandpaper to take these settings.

Use a melting pot for each alloy

Check frequently the state of the melting pot; small cracks are normal.

When the cracks become crunches throw the piece away.

Do not leave casting residuals inside the melting pot.

Do not let the metal solidify inside the melting pot; always use the receptacle to use the residuals in a better way.

In the noble metal casting there is the necessity of using graffiti melting pots.

Stuck

When there is use of metal residuals it is common to have stuck of the metal in the melting pot.

The metal that is near the bottom casts while the pieces that are above are trapped to the walls of the melting pot. Interrupt the casting and with a spatula distribute the metal in a more homogeneous way.

10 Maintenance

Turn off the general switch before any maintenance.

- Every 6 months disconnect the water entrance hose and clean the sieve.
- Clean the equipment with a wet cloth.
- Clean the centrifuge chamber regularly removing the casting residuals and revetment.
- Check regularly if there is not any metal residuals stuck to the coil, especially after casting in which there is metal asperse.
- Clean the view finder sieves periodically.

Problems	Causes	Solutions
On the display it appears the message "WATER FAILURE"	Low pressure in the line The exit hose is strangulated	Check the pressure of the line. Check if there is not other opened tap in the same line. Check if the exit hose is not strangulated.
On the display it appears the message "OVERHEATING FAILURE"	The room temperature or water refrigeration is too high. The ventilation orifices in the base are blocked.	Provide a better ventilation in the room. Unblock the ventilation.
On the display it appears the message "OPENED DOOR FAILURE"	The door was not close properly	Close the door correctly.

Problems	Causes	Solutions
On the display it appears the message "COIL FAILURE"	Coil lowered in the beginning of the cast. Coil lowered in the preheating Coil lifted during the centrifuge	Correct the operation error
On the display it appears the message "TIME FAILURE"	The maximum time of 2 minutes for casting was exceeded.	This function is useful to avoid forgetting and the overheating of the metal.
On the display it appears the message "CENTRIFUGE ERROR"	Too high supply tension. Defect in the centrifuge system.	Check the supply tension. If the defect remains contact the technical assistance.
On the display it appears the message "Erro 1"	Low current in the power circuit	Check the network tension. Check if there are metal residuals in the coil. Check if the coil had any damage.
On the display it appears the message "Erro 2"	Overcurrent in the power circuit.	Check the network tension. Check if there are metal residuals in the coil. Check if the coil had any damage
The metal does not heat or takes too long to cast.	The power key is in the minimal position. Low tension Excess of metal	Rotate the power switch totally to the right. Check the weight of the metal.
There is water leakage on the workbench.	The water exit hose is strangulated.	Check the exit hose.

Problems	Causes	Solutions
Although the coil gauntlet is lowered the centrifuge does not start.	The gauntlet is blocked before the releasing cam inferior position.	Force the gauntlet down and check what is blocking it.
There is metal leakage on the melting pot.	The melting pot orifice is blocked. The acceleration is too high for the volume of metal used. There is excess of metal.	Unblock the orifice Reduce the acceleration and the rotation. Check the weight of metal that is used.
There are constant failures in the casted pieces.	Low acceleration Low rotation Low temperature	Correct the parameters

Technical information:

Supply Tension:-----220 volts bifásico
Network frequency:-----50/60 Hz
Maximum power:-----3.000 Watts
Maximum current:-----14 amp
Maximum casting capacity:-----80g
Casting maximum temperature:-----1.750°C
Water flow at 2kg/cm²:-----4 litros /min
Argon flow :-----10 litros / min
Maximum refrigeration
water temperature:-----35°C
Maximum room temperature.-----40°C
Net weight:-----63 Kg
Width:-----0,6m
Height:-----0,53m
Depth:-----0,72m

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