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Dear customer,

The EDG Equipamentos has more than 100 Authorized Technical Service Center spread over Brazil so that you can enjoy their service promptly and calmly. Besides offering all the services by means of the Authorized Technical Service Center, the EDG Equipamentos also provides a range of other benefits like courses for the use of the equipment and technical and scientific advisory body, always aiming for the customer satisfaction.

With authorized and trained professionals from the Factory, tools and proper equipment, the EDG Authorized Technical Service Center is characterized by providing a fast and perfect Technical Service Center in order to please the customer's expectations.

To our customers that are not in Brazil, in case they need, we ask to contact the closest International Representative.

We remind you that the correct use of the equipment associated to maintenance and care periodically taken contribute to the maintenance of your equipment, considered, moreover, as an indispensable condition for the warranty right.

We thank for you trust and preference by our products and we feel honored to provide one more EDG equipment. Now you are also part of our family that doesn't stop growing.

Sincerely,
EDG Equipamentos

INSTALLATION

INSTALL THE FURNACE IN A WELL-VENTILATED PLACE, ON A BENCH-TYPE TABLE, MOVING AWAY CURTAINS AND INFLAMMABLE MATERIALS.

Verify if the furnace specified tension is compatible with the power grid tension (110 volts). Verify if the electric wire gauged is compatible with the furnace's power. Follows in the table below the minimum gauge for each type of furnace and the recommended fuse.

Model	220 Volts
3000	Fuse 10 A Wire


Use an individual outlet to plug in the equipment; never utilizes any kind of multiple inlet. If possible, install the fuse before the outlet that will feed the equipment.

SPECIAL CARE

Special Care should be taken with flasks, resistance and door

- FLASK** Manufactured with ceramic material, thermal shock must be prevented by following the instructions below:
- Never try to accelerate cooling leaving the door opened.
 - When opening the furnace's door, while it's still hot, keep it opened as short time as possible in order to prevent cracks.
- RESISTANCE** To extend the resistance life, observe the following instructions:
- Prevent temperatures above 1.100 °C
 - Never place inside flasks materials containing tin, lead or chemical products
- DOOR** Cares that must be taken
- Always close it gently
 - Periodically apply graphite powder on hinges.
 - To open, pull the lever upwards and then open the door.
 - In order to close the door, lean it carefully against frame and push lever downwards.

CONTROL 3P-S

- 1- DISPLAY
- 2- INCREASE VALUE KEY
- 3- DIGITS CHANGE KEY
- 4- PROGRAM-ON LED
- 5- PROCESS-ON LED
- 6- ST ART / STOP (S/S) KEY
- 7- MAIN SWITCH TURN ON/ TURN OFF
- 8- PROGRAMMED TEMPERATURE LEVEL LED
- 9- TEMPERATURE PROGRAMMING KEY
- 10- PROGRAMMED VELOCITY LEVEL LED
- 11 - VELOCITY PROGRAMMING KEY
- 12- PROGRAMMED STEADY PHASE LEVEL KEY
- 13- STEADY PHASE PROGRAMMING KEY

OPERATION

Note:

"Level" is the time in which the furnace will remain at the programmed temperature.

EXAMPLES OF 3 LEVEL PROCESS

To a better illustrate, we will give a programming example to the following burning process:

Level 1 - From room temperature, heat up to 300°C at a heating velocity of 15°C/min, remaining at this temperature for 30 minutes.

Level 2 - From the 300°C, heat up to 950°C at a heating velocity of 20°C/min, remaining at this temperature for 10 minutes.

Level 3 - From the 950°C, it cools down to 850°C in a controlled velocity of a 5°C/min cooling, remaining at this temperature for 15 minutes.

	Temperature (°C)	Velocity (°C/min)	Steady Phase (min)
Level 1	T1 = 300	V1 = 15	P1 = 30
Level 2	T2 = 950	V2 = 20	P2 = 10
Level 3	T3 = 850	V3 = 5	P3 = 15

Note:

The presented parameters are illustrative, not corresponding to any specific metal.


Temperature: Programmed Temperature

Velocity: Heating or Cooling velocity

Steady Phase: Steady Phase Level time

Initially, the system will show the furnace internal temperature.

For the example we will adopt 25°C

ABOUT THE KEYBOARD

1 - To reset or type a number - Resetting: Press continuously the key \uparrow (2) until the number 0 blinks on the display. Press \uparrow key until you get the desired value.

2 - To the next digit: Press the button \rightarrow (3) until the display starts to blink in the desired position. In case of desired digit is leapt, press continuously the bottom \rightarrow until you get back to the desired position.

To turn back to the first digit, you should process the same way as described above.

3 -If the number blinking is the desired value you don't need to program it again. Press the button " \rightarrow " to pass to the next digit.

PROGRAMMING

Sequence:

- A - Temperature 1, Temperature 2, Temperature 3
- B - Velocity 1, Velocity 2, Velocity 3
- C - Steady Phase 1, Steady Phase 2, Steady Phase 3.

1. Turn on the Main Switch Key (7)
 2. The Display (1) will show the real furnace temperature (our example 25°C) and the LED "PROGRAM" (4) will lit indicating that system will be ready to receive the desired programming.
 3. Press "T" key (9) until the "T1 " Led is lit. The first display digit (thousands digit) will blink
 4. Press the key "↑" (2) until the blinking number is "0" (zero)
 5. Press the → key (3) once, so blinking digit moves to the hundreds column.
 6. Press the key "↑" until the blinking number is "3" (three).
 7. Press the → key (3) once, so blinking digit moves to the tens column.
 8. Press the key "↑" until the blinking number is "0" (zero)
 9. Press the → key (3) once, so blinking digit moves to the unit column.
 10. Press the key "↑" until the blinking number is "0" (zero)
 11. Press the key "T" once, until the "T2" Led is lit. Repeat the steps from 4 to 10 above, but substituting the values to: "0"; "9"; "5"; and "0".
 12. Press the key "T" until the "T3" Led is lit. Repeat the steps from 4 to 10 above, but substituting the values to: "0"; "8"; "5"; and "0".
 13. Press the key "T".
- The 3 levels temperature is programmed.

Note: The maximum temperature accepted by the control is 1200 °C (others under consult). If a higher temperature was programmed, the equipment turns automatically to 1200°C.

14. Press the key "V", until the "V1" Led is lit.

The second digit on the display (hundreds digit) will blink.

Repeat the steps from 6 to 10 but with the "V" button and values ("0", "1" and "5") for the first velocity, ("0", "2" and "0") for the second velocity and ("0", "0" and "5") for the third velocity

15. Press the button "V"

The three level velocities are programmed.

Note: The maximum heating velocity that the control accepts is 30°C/min (others under consult). If a higher velocity was programmed, the equipment turns automatically to 30°C/min.

The maximum cooling velocity is the normal furnace heating loss; therefore, there isn't any device to accelerate cooling.

PROGRAMMING THE STEADY PHASE TIME IN EACH LEVEL

16. Press the Key "P" (13) until the "P1" Led is lit

The first display digit (thousands digit) will blink.

Repeat the steps from 4 to 10 but with the "P" key and values ("0", "0", "3" and "0"), ("0", "0", "1" and "0") and ("0", "0", "1" and "5").

17. Press the "P" key.

The three steady phase level is programmed.

Note: The permitted maximum value is 9,999 minutes (166,5 hours or 6 days and 22,5 hours) for each level.

At this stage the furnace will be programmed and ready to operate.

To initiate process, press the key S/S (6).

The "PROCESS" Led will lit, indicating that this process is working.

The "T1" and "V1" LED's will lit indicating the furnace starts to heat in order to reach "T1" temperature and "V1" velocity programmed.

When it happens the "P1" Led will lit indicating that equipment is keeping the temperature during the programmed time to the level 1.

At the end of the steady phase time in level 1, the equipment will move automatically to the second level and so on, until the end of process when a beep will go off and the word "END" will appear on the display.

Note: In the prosthesis furnaces, when finishing the process the equipment will keep the previous temperature.

In order to cool it off, press the button S/S or turn off the equipment.

For other purposes the equipment will start to be cooling when finishing the process.

OBTAINING MORE OF THE EQUIPM

1- REUSING PREVIOUS PROGRAMMING

The data of the previous programming are kept in its memory. To retrieve it push S/S key.

Note: The program is memorized only if performed at least once.

2- VERIFYING PROGRAMMING ON MEMORY

Press desired parameter key (T, V or P)

Example: What is the programmed heating velocity at level 2?

Press key V twice until indicating V2 LED is lit. The number shown on display is the programmed velocity at this particular level.

In order to return to the burning process, press the parameter key verified until the PROCESS LED (5) is lit.

3- TIMER UTILIZATION

This function is very useful when there is the need of performing the heating process without an operator. It allows that the furnace can be programmed in advance and begin to operate effectively from the desired time.

In order to use the timer function, it is necessary that the programming is made where the Steady Phase temperature level 1 is slightly above the room temperature. The Steady Phase time should be programmed at level 1 for determine the number of minutes the system will automatically enter in heating mode.

This time can go up to 9,999 minutes. (166.5 hours or 6 days and 22.5 hours)

	TEMPERATURE (°C)	VELOCITY (°c/min)	STEADY PHASE (min)
LEVEL 1	T1 = 30	V1 = 5	P1 = 360
LEVEL 2	T2 = 1150	V2 = 10	P2 = 120
LEVEL 3	T3 = 100	V3 = 5	P3 = 45

Note 1 : Level 1 temperature higher than room temperature

Note 2: In practice, the system will operate with 1 or 2 temperature levels, being the first level reserved to the timer.

- Once programming is completed the process is started by pressing key S/S
- The system will execute the cycle 1.
- As T1 temperature is rapidly reached (30%), the system immediately starts the steady phase, initiating the time counting.
- After 6 hours, (360 minutes for a programming time P1), cycles 2 and 3 will be executed respectively.

4- BLOCKED PROGRAMMING

In process where is needed just 1 or 2 temperature levels, its possible to program that blocks non-desired levels.

This resource facilitates programming in such cases. For this reason, its important to initiate the programming ALWAYS FROM THE CYCLE 1 and programming the parameters level that should not be executed with value 0.

Example:

T1 , V1 P1 = 0 Do not execute programming

T2, V2, P2 = 0 Executes level 1.; cancel 2 and 3

T3, V3, T3 = 0 Executes level 1 and 2; cancel 3

5- STEADY PHASE MONITORING

In order to verify the remaining time to the level that has been executed, press the button " " (3). The remaining time will be shown for 2 seconds after that the display returns to show the temperature automatically.


6- PROGRAM CHANGING DURING OPERATION

In any given moment during the heating process the programming can be modified, by pressing the corresponding parameter key. The LED will lit at the desired level; modify the value with the keys → and ↑ .

In order to return to the burning process, press the modified parameter key till the PROCESS LED lit.

Note: This modification will not be kept on memory.

SPECIFICATIONS

MODEL	Dimensions						Volume Liters	Net weight Kg
	Outside (mm)			Inside (mm)				
	A	B	C	D	E	F		
3000	300	400	470	150	100	200	3.00	16.7

MODEL	Electrical data	
	Output Watt	Voltage volt 50/60 Hz
3000	1700	220

FEATURES

- Low mass. Allows high thermal efficiency.
- High heating and cooling rates.
- Compact.
- Modular construction. Easier maintenance.
- Low cost

Temperature controls

3P-Manual

- Set Point
- Alarm bell.
- Adjustable power with 3 levels.
- Heating pitot lamp.
- Protection against sensor breakage.
- Accuracy class 1.5%.
- Protection door system.

Applications

- Dental prostheses.
- Ash and charge assaying.
- Burning.
- Heat treatment.
- Assaying of suspended solids.
- Ignition tests.
- Miscellaneous analysis.

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